


Baby Respiration Sensor



- ◆ For use with the Sims Graseby MR10, MR10+ or MR10s Respiration Monitor and Infant Flow Advance

- ◆ Accurately monitors respiration rate.

- ◆ The non-kink plastic tubing is non-electrical.

- ◆ Performance is not impaired when covered by the baby's clothes and does not restrict movement.

- ◆ Latex free. Tube Length 100cm.

- ◆ Individually packed in cartons of 50.

- ◆ Single use.

Vio Healthcare

Vio Healthcare Ltd. BN27 2BY - England
Tel: +44 (0)1323 446130
Fax: +44 (0)1323 847417
Email: sales@viohealthcare.com
Website: www.viohealthcare.com