


Loopies


Double Vessel Loop

Loopies enable retraction, identification and occlusion of vessels during surgical procedures.



'Highly versatile dual mode feature gives choice of sizes in each sterile pack.'

Mode 1



Mode 2



- ◆ Individually sterile packed in cartons of 50
- ◆ Elliptical cross section
- ◆ Medical grade silicone
- ◆ 400mm Length
- ◆ Radio-opaque

Catalogue Numbers

	RED	BLUE
Box with 50 double minis	VX120r	VX121b
Box with 50 double maxis	VX130r	VX131b

Vio Healthcare

Vio Healthcare Ltd. BN27 2BY - England
Tel: +44 (0)1323 446130
Fax: +44 (0)1323 847417
Email: sales@viohealthcare.com
Website: www.viohealthcare.com