

Surgical Skin Markers

- ◆ Specially formulated NON-TOXIC skin dye.
- ◆ Marks remain visible for several days.
- ◆ Choice of two line widths.
- ◆ Gamma Sterilized in fibre-free peel pouches.
- ◆ Big barrel for better grip.
- ◆ Individually sterile packed in cartons of 50.
- ◆ Private labelling for large volumes.
- ◆ Single use.



Vio Healthcare

Vio Healthcare Ltd. BN27 2BY – England
Tel: +44 (0) 1323 446130
Fax: +44 (0) 1323 847417
Email: sales@viohealthcare.com
Website: www.viohealthcare.com

The ideal Skin Marker for withstanding alcoholic skin preparations

As confirmed by tests in an article in British Journal of Plastic Surgery
(Volume 54 Number 6 2001)

**Allows the surgeon to demarcate selected sites,
prior to or during surgical procedures**

Non-toxic, violet ink is scrub proof

**The skin is usually marked and left to dry for at least
60 seconds prior to scrubbing**

**REF: VX100 Bold Line
REF: VX110 Bold Line + Ruler**

The marking instrument for

- ◆ Varicose Veins
- ◆ Thyroidectomy
- ◆ Orthopaedics
- ◆ General Surgery
- ◆ Radiology etc

**REF: VX200 Fine Line
REF: VX210 Fine Line + Ruler**

The marking instrument for

- ◆ Plastic Surgery
- ◆ Micro Surgery
- ◆ Intricate general surgery
facial, hand, etc
- ◆ Paediatric surgery